Figure 17A

Yield /g tissue Sterility (Bactec)

Attorney Docket No. 069952-0201 **In-Process Test** Comments Concentration in RC-3B 800g centrifuge 5°C Collect pellet and resuspend in RPMI Enrich for viable cells by Layer 100 mL RPMI mixing with 25% Iodixanol (1:1), and centrifuging in a Centrifuge @ 2000 rpm Cobe 2991 15 min **Pooled Mixing Band** Viability Collect the band at the top and pool Cell count & conc. Post washing Pellet Wash the Band twice with RC-3B centrifuge, 800g for RPMI 1640 (w/o phenol red) Viability and collect pellet 10 min, 5°C Re-suspend in 1/2 the volume of HTS required to achieve a nominal cell Viability concentration of 6 x 10⁷ cells/mL Add remaining volume of HTS Viability required to achieve nominal cell concentration (6 x 10^7 cells/mL) Cell count & conc. Fill cryobags with 1.5 mL of cell suspension Select a sample for: Add 1.5 mL of cryobuffer (concentrated) Viability Freeze using a programmable freezer (Cryomed) Cell count & conc. Store in Vapor phase liquid nitrogen freezer, Sample for release testing after 24 h Sterility: 3 bags from beginning, middle and end of filling

Figure 17B

1 bag at random

I bag at random

Endotoxins:

Characterization: